



60 PL110 00

PRESSLAM 1.1

Automatic mounting press - simple or double mounting

Description

- Automatic mounting press, new generation.
- Easy and friendly to use.
- Fully automatic, electronic controlled.
- All heating parts of the machine are covered and secured.
- Opening/closing system of the chamber by manual secure half turn lever.
- Pre-heating function to decrease the mounting time.
- Control and posting of parameters by color touch panel.
- 300 programs to save with free configuration.
- Saving/loading of programmes by USB and Ethernet connection.
- Select 4 different pressure methods.
- Automatic water cooling (with water saving function).
- Interchangeable heating and cooling unit (not included).
- Easily interchangeable mold assemblies (not included)
- Capable of double mounting.
- Solid steel case construction, power coated.





Technical data

- Mould assembly	ø 25 to 50mm (not included)
- Heating power	4x600W
- Temperature range	20 - 200°C
- Pre-heating fonction	40°C (disable without using of the machine)
- Heating cycle	0 - 30 min
- Pneumatique pressure	0 - 1600 daN
- Cooling cycle	0 - 30 min
- Air supply	6 to 10 bars max.Filtered 10µm
- Electrical supply	230 Volts single phase - 50 Hz

Dimensions & weight

250x480x620mm - 54 kg

Included

60PL110222 User and maintenance guide for presslam 1.1 in english

Option (has to be ordered with machine)

600061190 Chinese power plug

Code

Additional standard materials

60 PL2H1 00 Heating and cooling unit for mould assembly **ø25, 30 and 40mm (1-1/4" & 1-1/2")**

60 PL2H2 00 Heating and cooling unit for mould assembly **ø 50mm or 2"**

Molds assembly, single and double mounting.

60 PL2C1 00 mold assembly ø25,4mm

60 PL1C2 00 Mold assembly Ø30mm

60 PL1C4 00 Mold assembly Ø40mm

60 PL1C5 00 Mold assembly ø50mm

60 PL1C2 10 mold assembly 1-1/4"

60 PL1C4 10 mold assembly 1-1/2"

60 PL1C5 10 mold assembly 2"

60 PL2I1 00 Spacer ø25,4mm for double mounting

60 PL2I2 00 Spacer ø30mm for double mounting

60 PL2I4 00 Spacer ø40mm for double mounting

60 PL2I5 00 Spacer ø50mm for double mounting

60 PL2I2 10 Spacer 1-1/4" for double mounting

60 PL2I4 10 Spacer 1-1/2" for double mounting

60 PL2I5 10 Spacer 2" for double mounting

60 00375 90 Filtration kit for compressed air - 5µm